## 'Vortis' Insect Suction Sampler

Patent No. 9207468.1





Sole Manufacturers and Exporters

## Burkard Manufacturing Company Limited

UK and Export Sales Division Woodcock Hill Industrial Estate Rickmansworth Hertfordshire WD3 1PJ England

Telephone: 01923 773134/5 Fax: 01923 774790 E-mail: sales@burkard.co.uk Website: www.burkard.co.uk A lightweight portable suction sampling system designed to extract insects from herbage or debris with many distinct advantages over previously accepted methods.

- 1. The insects collected are immediately visible.
- 2. No nets, bags or filters.
- 3. Collecting vessel may be rapidly removed or replaced.
- 4. A wide range in size and shape of collecting vessel may be used.
- 5. Insects collected are safely deposited and may be returned to their habitat if required.
- 6. Light and fully portable system easily operated by one person.
- 7. NO losses from the catch when transferring collected material.

- 8. Ease and speed of collecting previously impossible to achieve.
- 9. Collected material does not impede air flow.
- 10. Air flow and therefore uniformity of collecting action can be monitored at all times.

This NEW design makes a radical break with the past in that the insects are deposited directly into a detachable transparent collecting vessel.

NO BAGS, NETS OR FILTERS ARE REQUIRED.

## **General Specification**

Air throughput	10.5m³/min.
Weight inc. Engine	7.8kg
Overall height	930mm
Collecting Area	0.2 sq/m